

CERTIFICATE OF ANALYSIS

Product	accu-beads®+	Part Number	710257
----------------	--------------	--------------------	--------

Concentration Type	Acceptable Range	Lot #		Date of Manufacture	Expiration Date
18 M/ml	± 2.5 M/ml	250210181		2025.10.02	2027.10.02
Test Type and Quantity	20 µm Fixed Coverslip Chamber Calibration	S.D.	CV%	Calculated Hemacytometer Calibration Range*	
1 L - Bio	18.20	0.28	1.56	23.28-24.02	
1 L - QC	17.79	0.68	3.82	22.24-24.01	
5 ml - QC	18.68	0.37	1.98	23.80-24.76	
QC Analysis by: Carlos A - 10/16/25		Reviewed by: Bryan Saint-Louis			
Date: 2025/10/16		Date: 2025/10/16			
Signature: Gary A. Wyor		Signature: Bryan Saint-Louis			

* Calculation based on hemacytometer flow dynamics ratio as described in accu-beads Research Report 45. Please contact Hamilton Thorne for more information.