

CERTIFICATE OF ANALYSIS

Product	accu-beads®+	Part Number	710258
----------------	--------------	--------------------	--------

Concentration Type	Acceptable Range	Lot #		Date of Manufacture	Expiration Date
35 M/ml	± 5.0 M/ml	250310351		2025.10.03	2027.10.03
Test Type and Quantity	20 µm Fixed Coverslip Chamber Calibration	S.D.	CV%	Calculated Hemacytometer Calibration Range*	
1 L - Bio	35.81	0.30	0.84	46.16-46.94	
1 L - QC	35.70	0.69	1.93	45.84-46.31	
5 ml - QC	37.46	0.48	1.29	48.06-49.32	
QC Analysis by: Carlos A. Ullas		Reviewed by: Bryan Saint-Louis			
Date: 2025/10/16		Date: 2025/10/16			
Signature: [Signature]		Signature: [Signature]			

* Calculation based on hemacytometer flow dynamics ratio as described in accu-beads Research Report 45. Please contact Hamilton Thorne for more information.