

UNONIUM³⁸® NEUTRALIS*

Disinfectant solutions adapted to your business

It is essential to use high-performance products which pass reproductive toxicity testing, in order to chart their total safety for gametes, and to reduce as much as possible any in-vitro stress which could alter embryo development.

UMONIUM38® NEUTRALIS is currently the only Medical Device-level range which meets the highest efficacy standards for all pathogens whilst also having passed the reference MEA (Mouse Embryo Assay) and HSSA (Human Sperm Survival Assay) tests.

The range can be used in all sectors, including within clinical units, without chemical risks or risks of interference; It hugely optimises the protocol's chances of success.



Distributed by

HAMILTON THORNE

A disinfectant that suits each of your needs!



UNONIUM³⁸ NEUTRALIS

Concentrated solution for dilution at 0.5% or 2.5% for application by soaking or damp wiping.

1L of concentrate makes up 200L of diluted solution (0.5%).

EXAMPLES OF APPLICATIONS (CLEANING AND DISINFECTION)

- Equipment (from operating room, cryoconservation, etc.)
- Suction pump
- Imaging/ultrasound equipment
- Instrumentation



UNONIUM³⁸ NEUTRALIS SPRAY

• Ready-to-use spray with double nozzle (fine spray or foam) for application by damp wiping.

• 500 ml can treat a surface area of 100 m².



UNONIUM³⁸ NEUTRALIS TISSUES

• Ready-to-use wipes for application by damp wiping.

• One wipe can treat a surface area of 3 m².

**Contact our dedicated team to ensure a seamless supply
of essential tools for your IVF lab!**

 **Distributed by**
HAMILTON THORNE



Visit our website!