

CERTIFICATE OF ANALYSIS

Product	accu-beads®+	Part Number	710257
----------------	--------------	--------------------	--------

Concentration Type	Acceptable Range	Lot #		Date of Manufacture	Expiration Date
18 M/ml	± 2.5 M/ml	241712181		2024.12.17	2026.12.17
Test Type and Quantity	20 µm Fixed Coverslip Chamber Calibration	S.D.	CV%	Calculated Hemacytometer Calibration Range*	
1 L - Bio	18.51	0.08	0.42	23.97-24.17	
1 L - QC	18.81	0.25	1.34	24.7-25.80	
5 ml - QC	18.35	0.62	3.39	23.05-24.67	
QC Analysis by: Gaby C. Miller		Reviewed by: Nicole Thibodeau			
Date: 2025-03-06		Date: 2025-03-10			
Signature: Gaby C. Miller		Signature: Nicole Th			

* Calculation based on hemacytometer flow dynamics ratio as described in accu-beads Research Report 45. Please contact Hamilton Thorne for more information.